

Standards of SMBG

عضو فدراسیون بین المللی دیابت

Sima Abbasi

Diabetes Educator

Gabric Diabetes Education Association

Renovating Diabetes Education

2023

SMBG roles

Treatment decisions

Critical event detection

CGM double check

Closed-loop system calibration



Accurate SMBG

- Holistic Understanding
- Factors Affecting Accuracy

01

Pre-analytical Factors

02

Analytical Factors

03

Post-analytical Factors



Holistic Understanding of Accuracy




Pre-Analytical Factors



Any factor that can affect the reliability of a test result occurring before the sample is analyzed.

Pre-Analytical Factors

- Lot Specific Coding
 - Test Strip Underfilling or Overfilling
 - Poor Surface Cleansing:
 - Fruit handling
 - Use of Specific Lotions
 - Inappropriate Strip Storage
 - Extraneous fluid at the testing
- Pseudo-hyperglycemia**
- Pseudo-hypoglycemia**
- 

Coding: A Source of Error



Coding determines the relationship between the electrical signal produced by the strip and the reported blood glucose.

IMPORTANT

Approximately 16% of EDUCATED patients miscode their meters.

16% Miscode Meters^{1,2}

16% of educated patients miscode their blood glucose meters



Incorrect coding may lead to measurement errors of $\pm 30\%$ or more.

Miscoding Causes Insulin Dose Errors

The calculated probability of some miscoded meters causing insulin dose errors of 2 units was as high as 50%⁴

	± 1	± 2	± 3	± 4	± 5
Miscoded meters 49.6%		50.5%	22.3%	1.4%	.04%
No coding meters 35.4%		1.4%	0%	0%	0%

Using Low dose algorithm

Insulin dose errors may occur⁴

No Coding or Auto coding meters avoid this error.

Under Filling Protection

Test performed

Results based on
TNO criteria

Accuracy (max. 15% deviation
from hexokinase-method)

60% passed


st strips results in

Reproducibility (max. CV 10%)


83% passed

Haematocrit dependency
(range 0.35 - 0.50 L/L)
(max CV 10% at maximum
values < 6.5 mmol/L or 1
for glucose values < 6.5 m

Not Filled



Filled



Haematocrit dependency
by manufacturers

Underfilling protection
(max 10% from result at minimal
volume or error mark)



systems had

underfilling protection

Under Filling Protection




20% of people have difficulties in filling test strips results in inaccurate result and painful pricking.



- 1) Minimal blood sample
- 2) Blood Volume Control (BVG) or Underfilling detection
- 3) Second chance Sampling: the possibility of re-sampling

External pressure may lead to unreliable readings



Different external pressures led to **$\geq 10\%$ differences** in glucose concentrations in 5%–13% of the participants.

On average, blood sugar readings were **lower** when people put pressure on the finger.

Hand washing has always been a problem!

With a sample of $0.3 \mu\text{l}$, **$1 \mu\text{g}$ of glucose (the weight of a dust particle)** will increase the blood glucose by **300 mg/dl** .



DO'S



DONT'S



Pre-Analytical Factors



URGENT

- ⚡ The impact of these factors has increased because of the smaller sample volume required by newer meters.
- 💡 Technological advances have addressed many of the pre-analytical sources of reduced accuracy.
- 💡 Comprehensive patient education remains important to address remaining pre-analytical issues.

International Diabetes Federation
Diabetes Education Association

Education is a MUST!



Patient demonstration of SMBG to the diabetes educator or health care provider is critical.

Table 3.

Tips for Successful Self-Monitoring of Blood Glucose Teaching

Use simple and specific steps at the patient's level of comprehension.

Be sure the patient can demonstrate the steps for SMBG.

Give your patient written recommendations for frequency and times of testing and desired results.

Observe SMBG procedure at follow-up visits.

Ask the patient to assess the relationship of SMBG with exercise, food, medications, and stress.

Specify which SMBG values are most problematic (especially low blood glucose) and discuss solutions with the patient.

Acknowledge the patient for goals achieved with SMBG.

Analytical Factors



Analytical Factors



All devices must include results of these tests in their labeling.

ISO 15197:2013 and FDA blood glucose meter accuracy standards

Table 7.1—Comparison of ISO 15197:2013 and FDA blood glucose meter accuracy standards

Setting	FDA (248,254)	ISO 15197:2013 (255)
Home use	95% within 15% for all BG in the usable BG range† 99% within 20% for all BG in the usable BG range†	95% within 15% for BG ≥100 mg/dL 95% within 15 mg/dL for BG <100 mg/dL
Hospital use	95% within 12% for BG ≥75 mg/dL 95% within 12 mg/dL for BG <75 mg/dL 98% within 15% for BG ≥75 mg/dL 98% within 15 mg/dL for BG <75 mg/dL	99% in A or B region of consensus error grid‡

BG, blood glucose; FDA, U.S. Food and Drug Administration; ISO, International Organization for Standardization. To convert mg/dL to mmol/L, see endmemo.com/medical/unitconvert/Glucose.php. †The range of blood glucose values for which the meter has been proven accurate and will provide readings (other than low, high, or error). ‡Values outside of the “clinically acceptable” A and B regions are considered “outlier” readings and may be dangerous to use for therapeutic decisions (256).


Post-Analytical Factors



Post-Analytical Factors




- Poor User Numeracy
- Misleading Unit Conversion
- Misleading Manual Log
- Software Processing Errors






SMBG: Component of the Education/Treatment Programme



عضو فدراسيون بين المللي ديابت



Holistic Approach

Take Home messages

- ▶ Inaccurate SMBG has both clinical and economical consequences.
- ▶ Technological advancements is necessary but not sufficient.
- ▶ Communicate with patients about accurate SMBG.
- ▶ Patient education is a MUST!